



BagCo Pty. Ltd.

ABN 48 097 819 291

Telephone: (03) 9464 0911

Facsimile: (03) 8677 1771

THE ALL NEW SPECIMEN BAG


TearZone®



Step 1

Technician places specimen into the front zipper pocket marked "Specimen Only" and zips the bag closed.

Patient paperwork is placed in the back pocket marked "Paperwork Only."



Step 2

When the laboratory technician receives the specimen, the bag is easily torn open at the TearZone no perfs, no blades!



Step 3

The specimen easily comes out of the bag protecting the lab technician, specimen, and paperwork!

TearZone® is a Trademark of BAGCO, Registration No. 2,338,575. Covered under patent number 6,149,302 and other patents pending.

TearZone® Specimen Bags are:

- ***Easy to load***
- ***Safer to dispense***
- ***Leakproof (Test Pressure 30KPa)***
- ***Packed in convenient dispenser packs***
- ***And, low cost !***

For further information or complementary lab trial samples please contact Bruce at Bagco – 0403 227 759